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USB CAMERA FOR MICROSCOPE MU-130Plus

Thank you for purchasing the USB Microscope Camera.
Please read this instruction manual thoroughly before use to insure proper operation and long service.

FEATURES

- 1.3M pixel USB Camera designed for easy use on a microscope. Simply connect to USB port for video image display.
- Powered by USB port; No separate power supply required.

NOTES


- This is a precision instrument – Handle With Care. Use caution to avoid shock or excessive forces on instrument when in use or transport.
- Do not use or store in direct sunlight or areas with high temperatures or high levels of dust. Keep in location where temperature and humidity are stable.
- Keep lens surface clean from dirt and fingerprints. To clean the lens, wipe gently with lens cleaning cloth or soft cotton cloth.
- When not in use, store in supplied case in location which has low moisture in order to prevent corrosion or internal mold growth.

SPECIFICATIONS

Model	MU-130Plus
Sensor	1/3 CMOS
Pixel Count	130M Pixels
Frame Rate	15fps (1280×1024) , 26fps (640×512)
Interface	USB 2.0
Power Supply	Supplied by USB port
Supported OS	Windows 7, 8.1, 10, 11
Min. recommended PC requirements	● CPU : Intel Dual Core $\geq$ 2.8GHz ● Main Memory : $\geq$ 2GB ● Hard Drive : 1GB Free space ● USB Port : USB 3.0 ● CD-ROM Drive
Weight (Body+Lens)	approx. 240g
Accessories	● Software CD ...1 ● Bore Adapters: ϕ23.2→ϕ30.0mm...1 ● ϕ23.2→ϕ30.5mm...1 ● USB Cable (approx. 2m) ...1

PART NAMES & ASSEMBLY

- ① Remove the eyepiece from the microscope.
- ② Remove the Lens Caps from the Camera Body and the Lens.
- ③ Follow the diagram below to connect.

- ④ Connect the MU-130Plus Camera to the PC USB Port using the USB Cable.
- ※ Keep the MU-130Plus Camera connected to the USB Port and continue with the additional installation steps.
- ⑤ When the installation procedure is complete, close ImageView and disconnect the MU-130Plus Camera USB Connector from the PC's USB Port.

TROUBLESHOOTING

- Q. "UCMOS01300KPA" does not appear on the Camera List.**
A. Camera is not connected. Please make sure USB Cables are properly connected.
- Q. Display image is dark or blurry.**
A. Insure there is sufficient lighting using microscope lamp or accessory light.
- Q. Driver software does not install. Viewer software does not work properly.**
A. Please confirm PC operating system is compatible; contact distributor.
- Q. Shortcut does not appear on the desktop.**
A. Shortcut will only be created for user with administrator privileges. If you want to create a shortcut for user without administrator privileges, go to the application installation directory to create the shortcut.
(Note) For default installation, the application is installed in the following directories:
C:\Program Files\MII\ImageView
C:\Program Files(x86)\MII\ImageView

UNINSTALL PROCEDURE

- In the Windows Control Panel, Select "Uninstall a program" and then select ImageView and uninstall.

FEATURE SUPPORT

Viewer Software (ImageView) provides supports for the following features:

- Still image capture and save
- Video capture and save

- ※ Other features are not covered.
- ※ Please note, support for additional features can not be provided even if requested.

USING THE VIEWER SOFTWARE (ImageView)

※Screen appearance will vary for different PC settings.

Image Preview

- ① Double click on ImageView icon to start the application.
- ② Select "UCMOS01300KPA" from the list and the Preview Screen will show an image.
- ③ Focus the image following the instructions for your microscope.
- ④ For some microscopes the image may display left/right reversed. Go to the "Flip" menu and check ☒ "Horizontal" option to mirror the image.


Still Image Capture and Save

- ① Click the "Snap" button to create a new image tab.
- ② Click on the "Video" tab to return to viewing the Preview Screen.
- ③ To save an image, double click on the tab for that image.
- ④ Enter the file name, type, and destination, then click the "Save" button.

Quick Functions

Quick Save

The Toolbar Quick Save icon will save images using settings specified in the Quick Save Menu.

Quick Record

When "Record" button is clicked, recording will start without a dialog box, and will use the setting specified here.


Image Quality Settings

White Balance

- Automatically adjusts the color tone using a white area of the image on the Preview Screen as a reference.
- ※ Color tone can also be adjusted using Color Temp / Tint sliders, but this is difficult to do accurately. We recommend using the White Balance button.
- When the White Balance panel is open, the Preview Image will display a red rectangle. This is the "white" reference area and the entire image will be adjusted based on this.
- Press the "White Balance" button to automatically adjust.
- ⚠ The area of the red reference rectangle may appear white on the Preview Screen even though the sample has a slight color tint. In this case the White Balance can not be set using the sample. We recommend removing the sample and setting the White Balance off the illumination light. (For epi-illumination, use a sheet of white paper.)

Power Frequency (Anti-flicker)

- This setting prevents a rolling bar or flicker on the image which can occur if the microscope camera frame rate is synchronized with the light illuminating the sample.
- Set according to the type power used for lighting the sample.
 - AC (50Hz) : Eastern Japan, Europe, etc.
 - AC (60Hz) : Western Japan, North America, etc.
 - DC : When using DC power.

Video Capture and Save

- ① Click the "Video" tab to display the Preview Screen.
- ② Click the "Record" button.
- ③ Select the video format and encoder.
- ④ Click "Next" and enter file name and directory to save the video.
- ⑤ Click "Finish" to start recording.
- ⑥ Click the "Record" button again to stop recording.

- If additional settings are required, click "Next" to specify "Bit Rate", "Time Limit", "Frame Rate", etc.

- Check ☒ "Time Limit" and specify time in seconds to end recording after the desired period.

INSTALLING ImageView SOFTWARE

Administrator privileges for computer are required for installation.

- ① Double click on the file "MIIImageViewSetup.exe" in the CD-ROM folder.
- ② When the installation guide screen appears, click on "Next" to continue.
- ③ Click on the "Browse" button to specify the directory for installation; when done click on "Install" to start the installation.
- ④ Click on "Finish" when the installation is complete.
- ⑤ The icon shown to the right will be placed on the desktop when the viewer software is installed.