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



STEREO ZOOM MICROSCOPE WITH LED ILLUMINATION

Model No. SMZDL-0790

Thank you for purchasing the Zoom Stereo Microscope with LED Illumination SMZDL-0790.

Please read this instruction manual thoroughly before use to ensure proper operation and long service.

PART IDENTIFICATION & SETUP

SPECIFICATIONS

Rated Voltage	AC90-240V 50/60Hz
Magnification	7~45× (for 10× eyepiece) 14~90× (for 20× eyepiece)
Objective Lens	0.7~4.5× (Zoom Ratio 1:6.4)
Eyepiece Lens	WF10× (FOV = 20mm) • WF20× (FOV = 10mm)
Working Distance	100mm
Eyepiece Angle	45°
Focus Adjustment	Movable Body
Stand	Rack and Pinion Type
Stage Plate	Φ95mm × 2pcs ①Acrylic Plate ②Black/White (Plastic Plate)
Eye Width Adj. Range	55~75mm
Diopter Adj. Range	±5 diopter
Illumination	Upper Light (Epi-Illuminator) White LED Lower Light (Transmitted) White LED
Base	L275 × W205 × H40mm
Weight	approx. 4.2kg

NOTES

- This is a precision instrument – Handle With Care. Use caution to avoid shock or excessive forces on instrument when in use or transport.
- Use and store in location with low humidity and stable temperatures. Avoid locations with dust, direct sunlight, excessive heat, or high vibration.
- Keep lens surface clean from dirt and fingerprints.
- Clean the lens gently with lens cloth or soft cotton cloth.
- Do not turn the left and right Zoom Controls in opposite directions; this will damage the instrument.
- Attach the included Upper Light according to the procedure below:

Upper Light (Epi-Illuminator) Installation

- ① Attach the illuminator with the left and right screws.
 - Adjust the mounting position and angle by sliding it.
- ② Connect the illuminator power cord to the top of the main unit.

SAMPLE PLACEMENT

- ① Place the sample in the center of the Stage Plate.
- ② If required, use Stage Clips to hold sample in place.

EYE CUP USE

- For use without eyeglasses: Hold the Diopter Adj. Ring to prevent it from rotating, turn the eyepiece, align the eye cup with the observer's face, and bring the eye close to observe.
- For use with eyeglasses: Remove Eye cups to use.

EYE WIDTH ADJUSTMENT

Hold both sleeves in each hand and look through the Eyepieces. Adjust the Sleeves to align the left and right field of view seen through the Eyepieces.

FOCUS & DIOPTR ADJUSTMENT

- ① Rotate the Zoom Control for highest magnification.
- ② Look through the right Eyepiece with the right eye and adjust the focus on the sample using the Focus Knob.
- ③ Rotate the Zoom Control for lowest magnification.
- ④ Using the Right Diopter adjustment, adjust the focus if required.
- ⑤ Now rotate the Zoom Control back to highest magnification and check the focus. Repeat steps ②~④ to achieve the best focus accuracy.
- ⑥ Rotate the Zoom Control for highest magnification. Look through the left Eyepiece with the left eye, and adjust focus using the Left Diopter Adjustment Ring.

FOCUS KNOB TENSION ADJUSTMENT

If the Focus Knob movement is too stiff or too loose, adjust the knob tension by holding either side of the Focus Knobs in place and turning the other side.

Turn it clockwise to add the knob tension, then eventually the head position is fixed. Turn it counterclockwise to reduce the knob tension.



FUSE REPLACEMENT

When the fuse blows, purchase a commercially available 250V 1A (20 mm × Φ5) cartridge fuse and replace it.

NSK Niigata Seiki Co., Ltd.

5-3-14, Tsukanome, Sanjo, Niigata, Japan, 955-0055
Tel. : +81-256-33-5522 Fax. : +81-256-33-5518
MAIL intl.sales@niigataseiki.co.jp
URL <https://www.niigataseiki.net/en/>

<Upper Light (Epi-Illuminator) Installation>

⚠ CAUTION

- Only use as directed. Improper use may cause accident or injury.
- When not in use store with protective cover, in location that has low moisture in order to prevent internal mold growth.