

Core-Shell HPLC / UHPLC Columns for Proteins and Peptides

Find newer methods in our Biozen Portfolio


See pp. 210-230

Ultra-High Resolution and Performance

Introducing Aeris, a specialized line of reversed phase core-shell HPLC / UHPLC columns, built exclusively for the ultra-high performance separation and analysis of proteins and peptides.


These columns can provide improved resolving power, selectivity, throughput, sensitivity, column lifetime, and method flexibility compared to other fully porous and core-shell columns typically used for bioseparations.

Aeris WIDEPORE



Large pore optimized for intact proteins and polypeptides

Aeris PEPTIDE



Small pore optimized for peptides and for peptide mapping

The precise architecture of Aeris core-shell particles provides dramatic leaps in performance in two important ways:

1 The thin, porous layer, or “shell”, decreases the diffusion path length, thus reducing the time it takes for biomolecules to adsorb/desorb into and out of the particle.

2 Expert manufacturing combined with tight packing specifications and high particle density reduces losses in efficiency and performance due to band broadening.

Aeris Core-Shell Particle

- High particle density helps create optimal bed structure which reduces band broadening effects of Eddy Diffusion
- Ultra-high performance on HPLC and UHPLC systems alike
- Reduced diffusion path improves efficiency



Fully Porous Particle

- Less homogenous bed structure leads to performance loss
- Ultra-high performance limited to sub-2µm particles on UHPLC systems
- Diffusion path limits efficiencies



The result is:

- **3.6µm core-shell particles** that can perform like sub-2µm columns on both HPLC and UHPLC systems at a fraction of the pressure
- **5µm core-shell particles** allow scale up to preparative dimensions
- **1.7µm and 2.6µm core-shell particles** that can provide higher peak capacities compared to fully porous sub-2µm columns on UHPLC systems



To see our entire BioSeparations column and accessory portfolio, visit: www.phenomenex.com/biopharm

Selecting the Optimal Aeris Column for Your Applications

Aeris core-shell columns are designed for the separation of complex protein and peptide mixtures. Chromatographers can easily narrow down the column(s) that has a high probability of success for their separation by selecting from a variety of phase, pore size, and particle size options.

Aeris PEPTIDE

Recommended for the separation of low molecular weight peptides and for peptide mapping.

- **XB-C18 chemistry best suited for resolving peptides**
- **1.7 μm , 2.6 μm , and 3.6 μm particles for method development flexibility between HPLC and UHPLC systems**
- **5 μm particle for peptide purification**
- **Small pore optimized for peptide diffusion**


For increased resolving power, use a longer column, preferably a 250mm (or 150mm for the Aeris 1.7 μm XB-C18). Due to the lower backpressure of Aeris 3.6 μm , one can easily run 250mm columns on both HPLC and UHPLC systems, AND one can couple multiple 250mm columns together and run them inline for even better results. For maximum UHPLC resolution, the 150mm length Aeris 1.7 μm or 250mm length Aeris 2.6 μm columns are excellent choices.

Aeris WIDEPORE

Recommended for the separation of intact proteins and polypeptides.

- **XB-C18, XB-C8, and C4 phases for alternate selectivities**
- **3.6 μm particle for system flexibility**
- **Large pore optimized for fast protein adsorption/desorption**

Because of the reduced hydrophobicity compared to fully porous 300Å columns, one should start gradients with reduced organic concentrations compared to other columns to improve peak shape of polar proteins and peptides. Shallower gradients compared to other fully porous columns may be appropriate.



Material Characteristics

Packing Material	Total Particle Size (μm)	Porous Shell (μm)	Core Size (μm)	pH Stability	Temp Stability $^{\circ}\text{C}$	Pressure Stability bar
Aeris WIDEPORE	3.6	0.2	3.2	1.5 - 9	90	600
Aeris PEPTIDE	1.7	0.22	1.25	1.5 - 9	90	1000
Aeris PEPTIDE	2.6	0.35	1.9	1.5 - 9	90	1000
Aeris PEPTIDE	3.6	0.5	2.6	1.5 - 9	90	600
Aeris PEPTIDE	5	0.6	3.8	1.5 - 9	90	600


Aeris™ Core-Shell LC Columns for Proteins & Peptides

Maximize Resolving Power with Unique Wide-Pore 3.6 µm Core-Shell Particle

3.6 µm core-shell technology combined with inert surface chemistries and tight packing specifications results in Aeris WIDEPORE columns delivering exceptional resolving power at significantly lower backpressures. Chromatographers now have the ability to

generate higher quality data than typically produced by columns packed with fully porous particles for every protein analysis – on HPLC or UHPLC systems.

Performance Equivalent to Sub-2 µm Particle at Low Backpressure




Ultra-High Resolving Power on UHPLC Systems with Aeris PEPTIDE 2.6 µm Columns of 250mm Length

The Aeris PEPTIDE 2.6 µm core-shell particle was designed with one purpose in mind: to enhance the separation and maximize the peak count of complex peptide maps on UHPLC systems. Because the 2.6 µm core-shell particle reduces backpressure on UHPLC

systems while maintaining similar efficiencies to sub-2 µm fully porous particles, longer columns can be used to further maximize the separation power while still being well within the backpressure constraints of the instrumentation.

Increase UHPLC Performance with Aeris PEPTIDE 2.6 µm Columns



Aeris™ Core-Shell LC Columns for Proteins & Peptides

Ordering Information

Aeris PEPTIDE 1.7 µm Minibore Columns (mm)				SecurityGuard™ ULTRA Cartridges*
Phase	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
XB-C18	00B-4506-AN	00D-4506-AN	00F-4506-AN	AJ0-8948

for 2.1 mm ID

Aeris PEPTIDE 2.6 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges*
Phase	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	3/pk
XB-C18	00B-4505-AN	00D-4505-AN	00F-4505-AN	00G-4505-AN	AJ0-8948

for 2.1 mm ID

Aeris PEPTIDE 2.6 µm MidBore™ and Analytical Columns (mm)				SecurityGuard ULTRA Cartridges*	
Phase	150 x 3.0	150 x 4.6	250 x 4.6	3/pk	3/pk
XB-C18	00F-4505-YQ	00F-4505-EQ	00G-4505-EQ	AJ0-8947	AJ0-8946

for 3.0 mm ID

for 4.6 mm ID

Aeris PEPTIDE 3.6 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges*
Phase	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	3/pk
XB-C18	00B-4507-AN	00D-4507-AN	00F-4507-AN	00G-4507-AN	AJ0-8948

for 2.1 mm ID

Aeris PEPTIDE 3.6µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges*
Phase	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
XB-C18	00B-4507-EQ	00D-4507-EQ	00F-4507-EQ	00G-4507-EQ	AJ0-8946

for 4.6 mm ID

Aeris PEPTIDE 5 µm Analytical Scout and Semi-Prep Columns (mm)					SecurityGuard ULTRA Cartridges*	SecurityGuard SemiPrep Cartridges**
Phase	150 x 4.6	250 x 4.6	150 x 10.0	250 x 10.0	3/pk	10 x 10 /3pk
XB-C18	00F-4632-E0	00G-4632-E0	00F-4632-N0	00G-4632-N0	AJ0-8946	AJ0-9317

for 4.6 mm ID

for 10 mm ID

Aeris PEPTIDE 5 µm Axia™ Packed Preparative Columns (mm)			SecurityGuard PREP Cartridges†
Phase	150 x 21.2	250 x 21.2	15 x 21.2 /ea
XB-C18	00F-4632-PQ-AX	00G-4632-PQ-AX	AJ0-9318

for 21.2 mm ID

Aeris WIDEPORE 3.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges*
Phases	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	3/pk
XB-C18	00B-4482-AN	00D-4482-AN	00F-4482-AN	00G-4482-AN	AJ0-8783
XB-C8	00B-4481-AN	00D-4481-AN	00F-4481-AN	00G-4481-AN	AJ0-8785
C4	00B-4486-AN	00D-4486-AN	00F-4486-AN	00G-4486-AN	AJ0-8899

for 2.1 mm ID

Aeris WIDEPORE 3.6 µm Analytical Columns (mm)				SecurityGuard ULTRA Cartridges*
Phases	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
XB-C18	00D-4482-EQ	00F-4482-EQ	00G-4482-EQ	AJ0-8769
XB-C8	00D-4481-EQ	00F-4481-EQ	00G-4481-EQ	AJ0-8771
C4	00D-4486-EQ	00F-4486-EQ	00G-4486-EQ	AJ0-8901

for 4.6 mm ID



SecurityGuard ULTRA
Holder with cartridge



Cartridge Holder

*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

**SemiPREP SecurityGuard Cartridges require holder, Part No.: [AJ0-9281](#)

†PREP SecurityGuard Cartridges require holder, Part No.: [AJ0-8223](#)



For HPLC Column Performance Check Standards, see pp. 424-425



Increase lab safety with HPLC / UHPLC solvent protection, see SecurityCAP™ products on pp. 417-418

For UHPLC system connections, see SecurityLINK™ UHPLC fingertight fitting system on pp. 336-337

For more about SecurityGuard ULTRA, see p. 335

For Core-Shell Performance Enhancement Kit, see p. 421