

## A Cost Effective Supported Liquid Extraction (SLE) Solution

Quickly and easily improve your liquid-liquid extractions by following a short, automatable two step extraction process. Packed with Diatomaceous Earth, Strata DE is a great alternative to traditional SLE products such as Biotage® ISOLUTE® SLE+, Thermo HyperSep™ SLE, and Agilent® Chem Elut® SLE.

### SLE Protocol

<b>Pre-treatment:</b>	Combine 100 µL of spiked urine, 15 µL Campbell Beta-Glucuronidase (part number: DR2102), 35 µL 100 mM Ammonium Acetate (pH 4), and 150 µL of 100 mM Ammonium Bicarbonate (pH 10).
<b>96-Well Plates:</b>	Strata DE 400 µL Biotage ISOLUTE SLE+ 400 µL
<b>Part No.:</b>	<a href="#">8E-S325-5GB</a> (Strata DE)
<b>Load:</b>	300 µL pre-treated urine sample onto plate (apply vacuum or positive pressure to pull/push sample into sorbent if necessary)
<b>Wait:</b>	6 minutes
<b>Elute:</b>	3x 600 µL Dichloromethane/IPA (95:5)
<b>Apply:</b>	Vacuum or apply positive pressure at 5-10" Hg for 10 seconds
<b>Dry:</b>	Sample under slow stream of Nitrogen at 30 °C
<b>Reconstitute:</b>	100 µL 0.1% Formic Acid/Methanol (4:1) with internal standard

### Recovery Values and % CVs: Strata DE vs. Biotage ISOLUTE SLE+

Analyte	Strata DE		Biotage ISOLUTE SLE+	
	% Recovery	%CV (n=8)	% Recovery	%CV (n=8)
6-MAM	98	9	88	16
Alprazolam	104	10	98	11
Benzoylcegonine	88	6	98	11
Buprenorphine	93	7	102	15
Codeine	99	12	93	9
Diazepam	107	7	104	6
Fentanyl	85	5	94	8
Hydrocodone	104	11	93	11
Hydromorphone	95	9	93	11
Lorazepam	94	8	98	8
Methamphetamine	92	16	102	8
Morphine	98	12	94	12
Norbuprenorphine	101	11	92	11
Nordiazepam	100	9	92	8
Norfentanyl	113	7	110	11
Oxycodone	97	5	93	11
PCP	90	7	98	6

## A Fast Extraction of 25-OH Vitamin D<sub>2</sub>/D<sub>3</sub> from Serum

Strata DE provides a simple extraction method with time and cost savings across all 3 QC levels.

### SLE Protocol

<b>Pre-treatment:</b>	Dilute 200 µL of human serum* with 100 µL of 5% Ammonium hydroxide (w/v), add 25 µL of 25-OH Vitamin-D <sub>3</sub> - <sup>2</sup> H <sub>6</sub> (1 µg/mL) and mix.
<b>96-Well Plate:</b>	Strata DE 400 µL
<b>Part No.:</b>	<a href="#">8E-S325-5GB</a>
<b>Load:</b>	Pre-treated sample and wait for 5 minutes
<b>Elute:</b>	Sample with 600 µL MTBE by gravity, wait for 5 minutes
<b>Repeat:</b>	Elution step twice by gravity, and after the final elution, apply 5-10 Hg vacuum to finish elution
<b>Dry:</b>	40 °C under N <sub>2</sub>
<b>Reconstitute:</b>	200 µL 0.1% Formic acid in Water/0.1% Formic acid in Methanol (30:70)

### Accuracy and Precision

	QCL	QCM	QCH
<b>Target Conc. (ng/mL)</b>	6	50	80
	25-OH-D <sub>2</sub>		
<b>Mean Conc. Found</b>	5.92	53.0	80.8
<b>STDV</b>	4.09	2.21	5.55
<b>CV%</b>	6.90	4.18	6.86
<b>Accuracy (%)</b>	98.7	106	101
<b>n</b>	6	6	6
	25-OH-D <sub>3</sub>		
<b>Mean Conc. Found (ng/mL)</b>	6.59	52.7	87.2
<b>STDV</b>	0.50	1.74	5.50
<b>CV%</b>	7.62	3.30	6.31
<b>Accuracy (%)</b>	110	105	109
<b>n</b>	6	6	6

\* Double Charcoal-stripped human serum was used to prepare all standards and QCs

Comparative separation may not be representative of all applications.

# Strata™ DE Diatomaceous Earth SLE

## Available for Large Volume Samples and High-throughput Cleanups

### Tubes

Ideal for large volume cleanups such as Food and Environmental applications.



#### Ordering Information

Strata DE (Diatomaceous Earth SLE Tubes)		
Strata DE Tube	2 mL Capacity, 12 cc	20 mL Capacity, 60 cc
Maximum Sample Volume (after dilution)	2 mL	17 mL
Recommended Elution Volume	2x 5 mL	3x 20 mL
Part No.	<a href="#">8B-S325-KDG</a>	<a href="#">8B-S325-VFF</a>
Unit	20/pk	16/pk

### 96-Well Plates

Ideal for smaller volume, high-throughput cleanups such as Bioanalytical samples.



#### Ordering Information

Strata DE (Diatomaceous Earth SLE) 96-Well Plates		
Strata DE 96-Well Plates	200 µL	400 µL
Maximum Sample Volume (after dilution)	200 µL	300 µL
Recommended Elution Volume	2x 600 µL	3x 600 µL
Part No.	<a href="#">8E-S325-FGB</a>	<a href="#">8E-S325-5GB</a>
Unit	2/pk	2/pk



For accessories that are compatible with Strata DE Supported Liquid Extraction (SLE) Products, see pp. 79-82



For more information on Strata DE, visit [www.phenomenex.com/stratade](http://www.phenomenex.com/stratade)

Recommended volumes are the expected loadability for most samples, however, it may be possible to load more than the stated capacity without breakthrough of the sample.

Comparative separations may not be representative of all applications.