



耐薬品表

Medium	セリベッター	セリベッタープロ
Acetaldehyde		+
Acetic acid (glacial), 100%		+
Acetic acid, 96%		+
Acetone		+
Acetonitrile		+
Acetophenone	+	
Acetylacetone	+	+
Acrylic acid		+
Acrylonitrile		+
Adipic acid	+	+
Agar (60°C)	+	
Allyl alcohol	+	+
Aluminium chloride	+	+
Amino acids	+	+
Ammonium chloride	+	+
Ammonium fluoride	+	+
Ammonium hydroxide, 30% (Ammonia)	+	+
Ammonium sulfate	+	+
Amyl alcohol (Pentanol)	+	+
n-Amyl acetate		+
Aniline		+
Barium chloride	+	+
Benzaldehyde		+
Benzyl alcohol		+
Benzylamine		+
Benzylchloride		+
Boric acid, 10%	+	+
Butanediol	+	+
1-Butanol		+
Butylamine		+
n-Butyl acetate		+
Calcium carbonate	+	+

Medium	セリベッター	セリベッタープロ
Calcium chloride	+	+
Calcium hydroxide	+	+
Calcium hypochlorite	+	+
Chloroacetaldehyde, 45%		+
Chloroacetic acid		+
Chromic acid, 50%		+
Copper sulfate	+	+
Cumene (Isopropyl benzene)		+
Diethylene glycol	+	+
Dimethyl sulfoxide (DMSO)		+
Dimethylaniline		+
Essential oil		+
Ethanol	+	+
Ethyl methyl ketone	+	+
Formaldehyde, 40%	+	+
Formamide	+	+
Formic acid, 100%		+
Glycerol	+	+
Glycol (Ethylene glycol)	+	+
Glycolic acid, 50%	+	+
Hexanoic acid	+	+
Hexanol		+
Hydriodic acid	+	+
Hydrobromic acid		+
Hydrochloric acid, 37%		+
Hydrogen peroxide, 35%	+	
Isoamyl alcohol		+
Isobutanol	+	+
Isopropanol (2-Propanol)	+	+
Lactic acid	+	+
Methanol	+	+
Methyl benzoate		+

Medium	セリベッター	セリベッタープロ
Methyl propyl ketone		+
Mineral oil (Engine oil)		+
Monochloroacetic acid		+
Nitric acid, 10%		+
Oxalic acid	+	+
Perchloric acid		+
Phenol		+
Phosphoric acid, 85%		+
Piperidine		+
Potassium chloride	+	+
Potassium dichromate	+	+
Potassium hydroxide	+	+
Potassium permanganate	+	+
Propionic acid	+	+
Propylene glycol (Propanediol)	+	+
Pyridine		+
Pyruvic acid	+	+
Salicylaldehyde		+
Salicylic acid	+	+
Silver acetate	+	+
Silver nitrate	+	+
Sodium acetate	+	+
Sodium chloride	+	+
Sodium dichromate	+	+
Sodium fluoride	+	+
Sodium hydroxide, 30%	+	+
Sodium hypochlorite	+	+
Sulfuric acid, 10%	+	+
Tartaric acid		+
Urea	+	+
Zinc chloride, 10%	+	+
Zinc sulfate, 10%	+	+

メモ：セリベッターとセリベッタープロはふっ化水素酸の使用はできません。

この耐薬品表は、実際にテストで得た結果を文献値より優先させています。リストに記載されている薬品以外でも多くの有機溶媒、無機塩溶液、洗浄剤組織培養用培地などが分注可能です。本リストに記載されていない薬品については関東化学(株)までお問合せ願います。

seripettor®と seripettor® pro は BRAND GMBH CO + KG の登録商標です。

当社による技術上の記載はユーザーの皆様へ情報を提供し、お役に立つとするものです。テストのための条件で得られた一般的な数値や結果を、実際の使用にあたって適応しようとする場合には、当社の責任では負いかねる様々な要素が影響することをご理解いただき、また、当社による技術上の記載を基にした特定の要求ができないことをご了承ください。実際の使用における有効性につきましては、ご使用される試料、作業環境に大きく左右されますので、お客様ご自身で試験して戴くようお願いいたします。

セリベッター / セリベッタープロの更に詳しい情報はブランド社のホームページ (www.brand.de) または関東化学までお問合せください。

仕様・価格は予告なしに変更する場合がございますのでご了承ください。

日本総代理店



関東化学株式会社

103-0023 東京都中央区日本橋本町3丁目11番5号 マルサンビル (03) 3667-6471
541-0048 大阪市中央区瓦町2丁目5番1号 (06) 6231-1672
812-0007 福岡市博多区東比恵2丁目22番3号 (092) 414-9361
HOME PAGE : <http://www.kanto.co.jp> E-Mail : instru-info@gms.kanto.co.jp

